

Can You See the Music?

PDQ 1 - Good Vibrations

Grades: 5-8

Time: 5 -15 minutes

Subject: Physics

Topics: Sound, Wave Properties



Overview

Catch a wave by listening to some sounds, seeing how they travel, and capturing them with both low and high-tech tools.

Background



Sound waves are created by vibrations at the source of the sound, such as a drum beat. These vibrations are transmitted through a medium, such as air, and can be felt physically. This is how a **microphone** works! As **sound waves** impact the **microphone**, a component in it moves. This movement is converted into electrical energy that can then be viewed, recorded, stored, or played back! In this PDQ, you will "see" **sound waves** causing grains of rice to move as well as see a visual display of the sound's intensity using Google Science Journal and databot™.

Objectives

Understand and Recognize:

- That sound is transmitted in waves.
- A **microphone** (sound sensor) vibrates from **sound waves** and converts this vibration to electrical energy.

What You'll Need

- Cereal bowl or open mouth container to stretch plastic wrap across
- Rice Grains -about 1 Tbsp
- Metal pan or pot and a metal spoon
- databot™ + Google Science Journal + Arduino IDE

Important Terms

Microphone: A **microphone**, sometimes referred to as a mike or mic, converts sound into an electrical signal.

Sound Wave: Sound is a vibration that travels in waves through a medium, such as air (or water, wood, etc.) These invisible waves have specific attributes such as frequency and amplitude. When a **sound wave** interacts with your eardrum, the vibrations of the wave are converted into a sound that you "hear" thanks to that amazing brain of yours!

Vibration: Vibration is an oscillating (back and forth) movement, like a vibrating reed in a clarinet. This vibration results in a **soundwave** that then travels through a medium, like the air.

decibel (db): Sound intensity is measured in units called **decibels**. The **decibel** scale is logarithmic, which means doubling the **decibel** units does not double the output, it can increase as much as 100 times! Normal conversation is about 60 **dB**, a soft whisper around 30 **dB**, and a lawn mower about 85 **dB**.

Prep (5 mins)

- Upload the Sound Intensity GSJ program to databot™, and place databot™ next to the bowl.
- Connect Google Science Journal to databot™ and open the databot™ sound sensor intensity function

PDQ 1 (10 mins)

- Stretch plastic wrap tightly across your bowl and place about 20 grains of rice on the tight surface. *(You are about to hit the pot with the spoon closely to the rice. What do you think will happen?)*
- Bang the pot or pan with your spoon next to the bowl and rice and observe what happens.
- Experiment by hitting the pot louder and softer. *Is it possible to hit the pot right next to the bowl, but the rice doesn't move?*
- Go further by maintaining a steady level of sound but moving away from the rice. *Does the distance seem to change the effect of the sound?*

databot™ has a sound sensor, a **microphone, built in. **Microphones** work by picking up sound vibrations and then converting the vibrations to electricity. You should see the sound level being displayed on GSJ.*

The Unit of measurement is decibels. Some examples of decibel levels: a pin drop = 10; rustling leaves = 20; babbling brook = 40; conversation = 60; alarm clock = 80; motorcycle = 100; rock band = 110; thunderclap = 120.

- Look at the grains of rice and the comparative sound intensity reading. *Is there a decibel level below which the rice does not move?*
- Try to maintain a sound level but move the source away from databot™. *Does the intensity stay the same regardless of distance or does it change?*

Great Work!

PDQ 2 coming up! Ready. Set. Go.

Another PDQ!

Educator Resources

Prep

- If you have not used databot™ and Google Science Journal before, go through the setup and use procedure.
- Upload and test the databot™ sound intensity program with Google Science Journal. Test your display connection.
- Study the background information and familiarize yourself with the learning objectives and terms for this activity.
- Setup and test the rice grain experiment and conduct the PDQ yourself before conducting it for your class.

Objectives

Understand and Recognize:

- That sound is transmitted in waves.
- A **microphone** (sound sensor) vibrates from **sound waves** and converts this vibration to electrical energy.

NGSS

- NGSS PS4.A Wave Properties

Misconceptions

- Sound is not dangerous. (Highlight the decibel levels at which damage to hearing can occur in your discussion).

Guiding Questions

- Why can databot™ pick up the sound even when the rice does not move?
- What do you think is the average decibel level in our classroom when we are working quietly.
- What about when we are having discussion like this one?

Additional Resources:

CDC: What Noises Cause Hearing Loss?

https://www.cdc.gov/nceh/hearing_loss/what_noises_cause_hearing_loss.html

Explain that Stuff - Sound

<https://www.explainthatstuff.com/sound.html>

Wikimedia Commons

<https://commons.wikimedia.org/wiki/File:Gnome-mime-sound-openclipart.svg>

Creative Commons License:

Can You See the Music? by Robert O. Grover & Team databot™ is licensed under a Creative Commons Attribution 4.0 International License. Permissions beyond the scope of this license may be available at databot.us.com/contact.



One little cube. Science on the move.

+1 208-451-2281
contact@databot.us.com
databot.us.com